

## EXPERIMENT NO. 5

### AIM :


To study the variation of thermo E.m.f of a thermo couple with temperature and determine its thermo - electric power.

### APPARATUS REQUIRED :

1. Thermo - Couple
2. Micro - voltmeter

### CONSTRUCTION OF THERMO - COUPLE :

The couple employed for the experiment is a Copper Constantan couple. A single piece of Constantan wire is welded at its two ends with two separate but similar copper wires so as to form two copper constantan junctions. These two junctions are separately kept in two beakers filled with water, thermometer and stirrers. One of the beakers will be heated while the other is to be kept at room temperature.



FIGURE



THEORY :

When two junctions of a thermocouple are kept at two different temperatures a thermo emf will be developed in the couple. Note the temperature of the hot junction from thermometer and the corresponding thermo emf from the micro voltmeter.

PROCEDURE :

Make connections as shown in fig.  
A graph is plotted with EMF measured along the Y-axis and temperature along X-axis.

OBSERVATION :

Temp. of the Cold  $j^n$  at the beginning of exp. =  $30^\circ\text{C}$   
 Temp. of the Cold  $j^n$  at the end of exp. =  $30^\circ\text{C}$   
 Average temp. of the Cold  $j^n$  =  $30^\circ\text{C}$

RESULT :

Plot the readings of difference of temperature of hot  $j^n$  and Cold  $j^n$  in the X-axis and thermo emf in the Y-axis. Draw a straight line passing through the maximum no. of points. Determine the slope of curve.

Teacher's Signature .....

S.No	(Temp. of hot junction) - (Temp. of cold junction) ( $^{\circ}\text{C}$ )	Thermo Emf (mv)		
		V(increasing)	V(decreasing)	$V_0$
1.	$30^{\circ}\text{C} - 30^{\circ}\text{C} = 0^{\circ}\text{C}$	0.133	0.130	0.132
2.	$35^{\circ}\text{C} - 30^{\circ}\text{C} = 5^{\circ}\text{C}$	0.142	0.140	0.141
3.	$40^{\circ}\text{C} - 30^{\circ}\text{C} = 10^{\circ}\text{C}$	0.150	0.156	0.153
4.	$45^{\circ}\text{C} - 30^{\circ}\text{C} = 15^{\circ}\text{C}$	0.167	0.175	0.171
5.	$50^{\circ}\text{C} - 30^{\circ}\text{C} = 20^{\circ}\text{C}$	0.192	0.192	0.192
6.	$55^{\circ}\text{C} - 30^{\circ}\text{C} = 25^{\circ}\text{C}$	0.200	0.208	0.204
7.	$60^{\circ}\text{C} - 30^{\circ}\text{C} = 30^{\circ}\text{C}$	0.219	0.224	0.221
8.	$65^{\circ}\text{C} - 30^{\circ}\text{C} = 35^{\circ}\text{C}$	0.235	0.244	0.240
9.	$70^{\circ}\text{C} - 30^{\circ}\text{C} = 40^{\circ}\text{C}$	0.260	0.255	0.258
10.	$75^{\circ}\text{C} - 30^{\circ}\text{C} = 45^{\circ}\text{C}$	0.281	0.270	0.276



Then the thermoelectric power = slope of this curve  
( $\tan \theta$ ) = thermoe.m.f developed /  $^{\circ}\text{C}$

For calculating the thermoelectric constant, two points  $(Q_1, C_1)$  and  $(Q_2, C_2)$  are chosen from the curve so obtained. Two quad. eq<sup>n</sup> formed are,

$$C_1 = aQ_1 + bQ_1^2$$

$$C_2 = aQ_2 + bQ_2^2$$

from these two eq<sup>n</sup> the values of  $a$  and  $b$  are determined.

Thermoelectric constant  $a = 4.89 \times 10^{-3}$  and  $b = 5.6 \times 10^{-5}$   
The neutral temperature  $t_n$  of the couple  $t_n = -a/b$ .

### PRECAUTIONS :

Water at the hot junction should be stirred constantly.

### CALCULATIONS :

$$C_1 = aQ_1 + bQ_1^2$$

$$\Rightarrow 0.24 = a(35) + b(35)^2$$

$$0.0685 = a + 35b \quad \text{--- (1)}$$

$$C_a = a\theta_a + b\theta_a^2$$
$$0.387 = 60a + b(60)^2$$

$$0.00545 = a + 60b \quad \text{--- (2)}$$

from (1) and (2)

$$\Rightarrow 25b = -0.0014$$
$$b = -5.6 \times 10^{-5}$$

Substituting  $b$  in eq<sup>n</sup> (1), we get

$$a = 4.89 \times 10^{-3}$$


The Neutral temperature  $t_n$  of the Couple,

$$t_n = -a/b$$

$$t_n = -\frac{4.89 \times 10^{-3}}{-5.6 \times 10^{-5}}$$

$$\Rightarrow \boxed{t_n = 87.32^\circ\text{C}}$$





SCALE

X-AXIS: 1 SSD = 1°C

Y-AXIS: 1 SSD = 20 mV

$$\text{SLOPE} = \frac{290 - 250}{50 - 37} = 3$$