Date	
0.0	

Expt. No	Page No.

	2
LXPERIMENT	

A	i	m	:
-	•		

To find the refractive index of liquid by travelling microscope.

Theory:

When viewed vertically from the air a point object inside the liquid; it appears to be raised by a small amount; depending on its depth below the surface and R.I. of the liquid relative to the air.

If the ray starting from (n) object appear to

If the ray starting from (p) object, appear to come from (p') which is the image of (p) object, appears to come relative to the air is given by H = u - op where op=(u) real depth and op'=(v) apparent depth i.e.,

H = Real depthApparent depth

Required Apparatus:

- 1> Travelling Microscope
- 2) Beaker (cross mark on the base)
- 3> Spirit level.
- 4) Tissue papper
- 5> Magnifying glass
- 6> Water

Teacher's Signature

1	Date		
	Page	e No	,,
	Real	Depth	Apparent

Expt. No.	
-----------	--

TABLE

Table for Experimental Data

All readings are in cm.

SiNo	Pan II	٦. د	2		•						
J. 140.	Main	ng of			s of 1			8 of	Q	Real Depth	Apparent
	Scale	Vernier 8 cale	Total	Main_ Scale	Vernier Scale	Total	Main	Vernier	Total	1 >	Depth
1.	5.1	1	5.101	5.8		5 000	Scale	Scale		(Q-P) u	(Q-P')V
		i×0.001		3.8	2 2x0.w1	5.802	7.8	2×0.001	7.802	2.701	2
		= 0.001			=0.002			= 0.002			
2.	5.1	1	5.101	5.9	3	5.903	8.3	4	8.304	3.203	2.401
		1×0.001			3x0.007		a change	4x0.001		0.203	2. 102
		= 0.001			= 0.003			= 0.004			
3.	5.1	1	5.101	6.0	2	6.002	8.6	6	8.606	3.505	2 /0/
		1×0⋅∞7			2×0.001			6×0.001		3.505	2.604
	-	= 0.001			= 0.002		1	= 0.006			
							1				
							P. Iv.				

Refractive Index of (given liquid) M = Real Depth
Apparent Depth


1) $u_1 = 2.701$ and $v_1 = 2$

Refractive index $H_1 = U_1 - 2.701 = 1.3505$

Teacher's Signature ___

t. No	Page No.
2 $u_2 = 3.203$	V ₂ = 2.401
then , RJ. (wakr)	$H_2 = U_2 - 3.203 - 1.3340$ $V_2 = 2.401$
3) u ₃ = 3.505	, V ₃ = 2.604
then, R.I. (waker)	$y_3 = u_3 - 3.505 - 1.3460$
Therefore mean R.I.	$Me = \underbrace{M_1 + M_2 + M_3}_{3}$
	= 1.3505 + 1.3340 + 1.3460 3
	= <u>4.0305</u> 3
R	$H_e = 1.3435$ of water = 1.3435

Teacher's Signature _



Date			

Expt. No	Page No.

	From	the	above	dasis	R.T.	of water	(Me)= 1.333		
1								1.343 + 1.333	= 1.338	
						2		2		

Result:

Hence the refractive of water calculated here is 1.338 (opprox.)

Percentage error = <u>Difference of true value Prakulate value x 100</u> true value


= 1.338 - 1.339 × 100 /.

1. error = 0.375 1.

Precautions: -

- 1> To take accurate reading use magnifying glass.
- 2) Do not alter distance of eye piece and object after taking initial reading.
- 3> After reading for error mark (p), do not adjust the adjustable screw of the eye piece, while reading of p' and q. Only adjust the adjustable screw of the vertical seale screw.

Teacher's Signature



Real Depth and apparent depth of a fixed point below the liquid